

Contour Recline E-Vertex

A world-first, fully electric procedure chair with touch control adjustable headrest for ease of operation, access and patient comfort.



Key Benefits

- **Full flexibility:** The infinitely adjustable headrest, range of chair heights with lie-flat and forward/reverse options places the surgeon and patient in exactly the right position
- **Safety and support:** Increased leg clearance, with an adjustable, rigid locking headrest and wrist support accessory, brings patients close to protect surgeons from neck and back strain.
- **Comfort and control:** Instantly increase patient comfort and enhance access with one-touch, electronic control of the scalloped or flat cushioned headrest; and pillow support accessory
- **Hygiene conscious:** Smooth, gap free platform aids cleaning and helps prevent fluid ingress to meet the highest infection control standards.

Application

Ideal for Eye, ENT, plastics and dental care.

Why Buy?

The Contour E-Vertex medical procedure chair has been designed with direct feedback and insights from the country's leading Eye, ENT, plastics and dental professionals. The breakthrough design combines state-of-the-art technology and ergonomics to deliver exceptional functionality and support for procedures requiring close proximity, control for the practitioner and patient comfort. With convenient, touch-button handsets and user-friendly mechanisms, allowing a surgeon and their support team to quickly, safely and smoothly adjust the headrest, chair height and position for the procedure, patient transport and recovery. [Also available as a Stretcher, view the Contour E-Vertex here.](#) *International patents pending

v-12112024

MODSEL

Making a world of difference.

modsel.com ►

Contour Recline E-Vertex



Low and safe

Excellent 460mm low height with fold-under drop-sides.

The ultra-low height is combined with drop-sides that glide down vertically under the chair top to provide safe access with no protrusions to catch on the patients legs or clothes.



Lies flat

The backrest and leg rest move independently or together from chair to bed.


You now have the flexibility of sitting a patient up to almost 90 degrees or lying down completely flat.



Smart transfer

Conveniently stored handles unfold for staff friendly transportation.

Lock in the double steer system, unfold the push handles and move forward effortlessly on the large castors.




Slim and clean

The ultra thin design safely hides all mechanical operation and has a one-piece, gap free top surface for easy cleaning.



Effortless control

Clearly marked headrest buttons allow for precise and supportive positioning of the patients head with ease.



Raise and lower

Adjust the height up and down while remaining level.

Specifications

- **Size:** 780 x 1150mm overall (1750mm lying flat)
- **Seat width:** 650mm with 100mm wider option available with 100mm thick premium mattress
- **Height Range:** 520 - 900mm with lie-flat and forward / reverse tilt
- **SWL:** 300kg
- **Castor Size:** 150mm twin wheel castors
- **Package Includes:** standard IV, vertical oxy holder, drop down footrest, meal tray, standard scalloped or flat headrest cushion option

Product Codes

- Contour Recline E-Vertex: MPC-C-REVE-E
- Contour Recline E-Vertex Package: MPC-C-REVE-E-P
- Pillow Support: MPT-CA-EPS
- Wrist Support: MPT-CA-EWS
- Flat Headrest Cushion Option: E0311

v-12112024

MODSEL

Making a world of difference.

modsel.com ▶