

# Dropside Options

There are a range of dropside options available to customise any stretcher to suit your special requirements.



## Key Benefits

- Multi-riser dropside options for closer spacings
- Split fold option for safer patient access
- Horizontal rail designs for longitudinal support
- Flush folding option for mattress transfer or x-ray
- Mid position lock option for easier patient access

v-07102022

# MODSEL

*Making a world of difference.*

[modsel.com](https://modsel.com) ►

# Dropside Options



## 16 riser dropside

Close up the spacing's to 100mm centers with extra risers for a safe side rail for patients requiring extra support.

Choice of folding from head end towards foot or foot end towards head (pictured).



## 7 riser dropside

Close up the spacing's to 250mm centers with extra risers for a safe side rail for patients requiring extra support.

Choice of folding from head end towards foot (pictured) or foot end towards head.



## Split folding dropside

Enhance access to either the head or foot end of the patient without compromising safety.

Fold all four sections down together or independently. Supplied standard with mid-position catches.



## Flush folding dropside

Fold down the dropside to flush with trolley top, allowing slide transfer of mattress or horizontal x-ray.

Uses the slim aluminium handrail section. Choice of folding from head end towards foot (pictured) or foot end towards head on any dropside riser configuration.



## Mid-position lock

Pair of catches with a mid locking position for special access to the patient.

Choice of folding from head end towards foot (pictured) or foot end towards head. Supplied standard on split dropside.



## Deluxe dropside

Full length glide-away siderails that provide horizontal support and a zero transfer gap.

## Why Buy

We understand that there are different ways that patients must be transported and procedures performed. Safety comes first, which is why we've developed a selection of side rail options, to ensure you can have the right outcome for your department's requirements.

## Product Codes

- 16 riser dropside - MPT-CA-16RD
- 7 riser dropside - MPT-CA-7RD
- Split fold dropside - MPT-CA-SDS
- Flush fold dropside - MPT-CA-FFD
- Mid position lock - MPT-CA-CMP
- Deluxe dropside - MPT-C-D

v-07102022

# MODSEL

*Making a world of difference.*

modsel.com ►