

End Rail

An end rail that can help to retain the mattress, patient or be used as a push handle.

Key Benefits

- Strong, 25mm diameter coated stainless steel rail
- 600mm wide x 200mm high
- Fits into foot end or optional head end accessory sockets

Product Codes

- End rail; MPT-CA-FER
- Head end accessory sockets; MPT-CA-ASH (pair)

Why Buy

Quickly drop this rail into the accessory sockets (optional accessory on the head end of a stretcher, excludes Contour Vertex) and use it as a higher strong point to move the stretcher. This rail has a variety of additional uses such as retaining a mattress or pillow from sliding off or a bar for a patient's feet to push against. These rails can be used on either end of the stretcher or both ends together.



Head end option



Both ends

v-20012023

MODSEL

Making a world of difference.

modsel.com ►