

Folding Monitor Board

A strong frame with push handles and folding board with an adjustable tie strap for retaining a monitor.

Key Benefits

- Powder coated steel frame
- Cross bars for patient feet protection
- One-piece hospital grade plastic board
- 50mm wide adjustable tie strap
- 450mm wide x 400mm deep tray surface
- 600 wide x 470mm high from mattress base
- 15kg weight capacity

Product Code

- MPT-CA-FMB

Why Buy

Add this removable monitor board accessory to the foot end of a stretcher to provide a solid push handle and fold-away board for securely mounting a monitor. The board hangs vertically out of the way and can be easily rotated outward and over to a stable, internal, horizontal rest position.



Fold over to a table



Writing surface option



Writing surface extended

v-20012023

MODSEL

Making a world of difference.

modsel.com ►