

Folding Monitor-Foot Board

A strong frame with both a folding foot and monitor board with an adjustable tie strap for retaining a monitor.

Key Benefits

- Powder coated steel frame
- One-piece hospital grade plastic board construction
- 600 wide x 410 deep x 540mm high overall size
- 520mm wide x 380mm deep monitor tray
- Monitor tray is 540mm high from mattress base
- 15kg weight capacity
- 50mm wide adjustable tie strap for monitors
- 500 wide x 250mm high foot board
- Adds 320mm to the length of the mattress

Product Code

- MPT-CA-FMFB

Why Buy

Add this removable monitor board accessory to the foot end of a stretcher to provide a solid push handle and fold-away board for securely mounting a monitor. The board hangs vertically out of the way and can be easily rotated outward and over to a stable, internal, horizontal rest position. Padded for patient comfort, the dropdown foot board folds with a simple 'lift and lay' action to quickly increase the length of the mattress area.



Fold over to a table



Fold down foot board



Writing surface option



Writing surface extended

v-20012023

MODSEL

Making a world of difference.

modsel.com ►