

# Dropside Padding

Provide extra protection for recovering patients with a selection of padded bumper systems for stretcher dropsides.



## Key Benefits

- Stiff board backed bumper system
- Lift and roll feature allowing use without removal
- Quick install pocket system
- 25mm high density padding
- Welded cover materials for easy cleaning
- Sizes to fit any dropside style

v-20012023

# MODSEL

*Making a world of difference.*

modsel.com ►

# Dropside Padding



## Padded Bumpers

Stiff plastic backed bumpers with welded stretch padding that resists high limb force.



## Unique access

Lift and roll to the outside of the dropside, allowing raise and lowering without removal.

Wall hanging point incorporated for easy storage.



## Padded Pockets

Tough vinyl padding on both sides of the dropside protecting both the patient and carer.



## Removable design

Unclip and quickly lift off for unrestricted patient access and easy cleaning.



## Sizing options

Both bumpers and pockets can be supplied as full length versions or to fit different dropside styles.

## Why Buy

A recovering or restless patient can cause themselves serious harm on the metal risers and rails of standard dropsides. We have developed some ways for the carer to quickly and easily increase the level of protection and provide a soft, padded surface for patients limbs. Our board system can be flipped to the outside of the stretcher allowing the dropsides to be raised and lowered without having to remove the bumpers from the trolley. Pockets can slide on over the dropside which are quick to use, however must be removed to lower the sides.

## Product Codes

- Padded Bumpers (pair); MPT-CA-PDB-P
- Padded Pockets (pair); MPT-CA-PDB-D

v-20012023

# MODSEL

*Making a world of difference.*

modsel.com ►