

# Contour E-Vertex

A fully electric stretcher with touch control adjustable headrest for ease of use, access and superior stability.



## Key Benefits

- **Full flexibility:** The infinitely adjustable headrest up and down tilt and excellent raise height allows the surgeon to place a patient in exactly the right position.
- **Safety and support:** The compact, rigid locking headrest with increased leg clearance brings patients close and keeps them stable while minimising surgeon neck and back strain.
- **Comfort and control:** Instantly increase patient comfort and enhance access with one-touch, electronic control of the scalloped or flat cushioned headrest; and lightweight removable frame with cushions to support a pillow during recovery.
- **Slim and clean:** The ultra-thin design hides all mechanical operation and has a one-piece, gap-free top surface for ease of cleaning and infection control.

## Application

Ideal for Eye, ENT, plastics and dental care.

## Why Buy?

We know precision care requires 100% confidence in the equipment being used. That's why the Contour E-Vertex combines state-of-the-art technology and ergonomic design to deliver exceptional support for the performance of procedures requiring close patient proximity, control for the practitioner and patient comfort. The slimline design of this extremely stable, rigid locking headrest provides excellent leg clearance – bringing the patient closer. While the touch button electronic positioning system\* allows a surgeon or their support team to quickly, safely and precisely position a patient's head. [Also available as a procedure chair, see the Contour Recline E-Vertex here.](#) \*International patents pending.

v-10072024

# MODSEL

*Making a world of difference.*

modsel.com ►

# Contour E-Vertex



## Slim and clean

The ultra thin design safely hides all mechanical operation and has a one-piece, gap free top surface for easy cleaning.



## Effortless control

Clearly marked headrest buttons allow for precise and supportive positioning of the patients head with ease.



## Raise and lower

Adjust the height up and down while remaining level.




## High tilt

Excellent raise height with up and down tilt adjustment.



## Low tilt

A practical low height with up and down tilt adjustment.



## Pillow Support option

Lightweight removable frame with cushions to help support a pillow during recovery



## Wrist Support option

Coated, 280mm wide, height adjustable and removable rest



## Fold-back design

The wrist support flips away to allow quick and easy access to the patient

## Specifications

- **Overall Size:** 2200 x 815mm (965mm wide option).
- **Patient Surface:** 2070 x 680mm (830mm wide option).
- **Height Range:** 515mm-885mm.
- **SWL:** 300kg.
- **Castor Size:** 200mm castors and Ezi-Trak steering.
- **Package Includes:** 75mm Premia Mattress, Vertical Oxy Holder (Rail Mount), Standard IV Pole, Foot Endrail, Pillow/Shoulder Support, 2nd Handset at Foot End of Stretcher & USB Port for charging of patient devices.

## Product Codes

- Contour E-Vertex: MPT-C-VE-E
- Contour E-Vertex Package: MPT-C-VE-E-P
- Pillow Support: MPT-CA-EPS
- Wrist Support: MPT-CA-EWS
- Flat Headrest Cushion Option: E0311

v-10072024

# MODSEL

*Making a world of difference.*

modsel.com ►