

# Dropdown Foot Board

Padded fold-down foot board to increase the length of the patient surface.

## Key Benefits

- High density foam covered in welded stretch material
- Tough board construction in a one-piece hospital grade plastic
- Stainless steel lift and lock folding hinge
- 530mm wide x 330mm high (250mm high cushion)
- Adds an extra 310mm to the length of a mattress

## Product Code

- MPT-CA-DDF

## Why Buy

Padded for patient comfort, the removable dropdown foot board has molded handles for pushing the stretcher and with a simple 'lift and lay' action, folds down to quickly increase the length of the mattress area.



Fold down

v-20012023

# MODSEL

*Making a world of difference.*

modsel.com ►