

Dropdown Push Handles

Strong fold-away handles with rubber grip for pushing the trolley.

Key Benefits

- Strong stainless steel construction
- Lift up and fold down
- Tough hygienic hand grips
- Tall head end handle option - 400mm high x 600mm wide (350mm between handles)
- Foot end handle option - 160mm high x 600mm wide (350mm between handles)

Product Codes

- Tall head end handles; MPT-CA-TDPH
- Foot end handles; PT-CA-DPH

Why Buy

To allow space for the many removable and common accessories used at the foot end of the stretcher, the handles are shorter to fold completely away with no obstruction to the accessory mounting positions. There is an option for a tall version of these handles to be mounted at the head end of the stretcher which provide high pushing points and ample clearance from the moving backrest.



Head end handles folded



Foot end handles



Foot end handles folded

v-20012023

MODSEL

Making a world of difference.

modsel.com ►